STERILE

EO

MD

**(UK) - INSTRUCTIONS FOR USE: - OESOPHAGEAL DOPPLER PROBES**

**1. Intended Use:**

The Oesophageal Doppler Probes are used with CardioQ™/CardioQ ODM™ Monitoring Systems. The Doppler Probes have a designated time limit and when that limit is exceeded the probes will cease to function. The I<sub>2</sub> Oesophageal Doppler Probes are for oral or nasal use in patients under general anaesthesia or sedated patients in ITU.

**Clinical benefit-** Proven effective for haemodynamic optimisation using 10% SVO to reduce complications and shorten length of hospital stay

**Indications-** The Deltex Medical CardioQ-ODM+ is indicated for use as a fluid management and cardiac output monitoring system

**Warning:** Possible expiration of probe whilst inserted into patient.

The user is constantly shown the usage time remaining and is notified when the probe has 20 minutes of usage time left, with both an audible and visual alert. This alert is again repeated when there is 5 minutes of usage time remaining. This should give ample time for the user to get a replacement probe if one is required.

Before use read the appropriate Operating Handbook.

Use only with the appropriate Monitoring System.

**2. Warnings - Do Not Use in the Following Situations - Patient Harm May Result:**

- Do not use after the expiry date printed on the shelf box and unit pack labels.
- Do not use undue force upon insertion. Remove if difficulty arises during insertion and seek advice.
- Verify coagulation status for the possibility of nasal bleeding prior to inserting.
- Use only by trained and qualified personnel, or under such supervision.
- Not recommended for use in patients under going intra-aortic balloon pumping.
- The probe may be unreliable during periods of aortic cross clamping.
- Any serious incident that occurs while using this product shall be reported to Deltex Medical Ltd and your member state competent authority.
- MR unsafe, do not use this equipment in the MRI scan room

**3. Cautions:**

Do not use if the pouch has been unintentionally opened, or if the probe or pouch show signs of damage. The probe shall be disposed of according to local procedures

Exercise caution when handling. The Probe contains an internal spring that will cause the probe to uncoil when released from the pouch.

Open the pouch at the chevron seal end.

Do not use scissors or sharp tools to open the pouch.

Do not use if the pouch has been opened, or if the probe or pouch show signs of damage.

**4. Contraindications:**

- Not for use on patients under 16 years of age
- Do not use with carcinoma of the pharynx; larynx; oesophagus.
- Do not use with aneurysms of the thoracic aorta.
- Do not use with proximal coarctation of the aorta.
- Do not use with tissue necrosis of the oesophagus or nasal passage.
- Do not insert nasally where there is nasal injury; where there is facial trauma; nasal polyps or where there may be the risk of brain injury
- Do not use in close proximity to laser surgery
- Not for use in patients with pharyngo-oesophago-gastric pathology and/or severe bleeding diatheses.

**5. General Information:**

The Doppler Probes are supplied STERILE.

Caution: Software compatibility: The Probe can only be used with an appropriate CardioQ/CardioQ ODM Monitoring System that has Release 5.60/1.03 or higher software installed. Probe output specifications is documented within the ODM Monitor handbook

**6. Removal of the Probe from the Unit Pack:**

Open the pouch at the chevron seal end and, holding the Probe and withdraw it from the pouch.

Allow the Probe to uncoil, keeping it free from contact with any object that may compromise its state of cleanliness.

**7. Insertion of the Probe into the Patient:**

Insertion should be carried out by an appropriately trained and qualified medical practitioner or by a person under the supervision of such.

Insertion of the Probe may be either orally or nasally. Note the Warnings above.

The Probe connector should be plugged into the Patient Interface Cable connector and the correct nomogram entered.

Apply a water-based lubricant to the tip of the Probe. This aids insertion and signal acquisition. Insert the Probe with the bevelled edge of the tip facing the back of the throat so that when in the oesophagus, the bevelled edge faces backwards and downwards towards the descending aorta.

Do not use excessive force upon insertion. Remove if difficulty arises and seek advice.

**Warning -** Risk of patient harm if forced.

The numbered depth markers on the Probe shaft are for use as a guide for oral and nasal placement. With oral placement, signal acquisition may normally be achieved at a depth of between 35cm (marker 1) and 40cm (marker 2) in an adult patient. With nasal placement, signal acquisition may normally be achieved at a depth of between 40cm (marker 2) and 45cm (marker 3) in an adult patient. The required depth of insertion will depend upon patient height/weight. Patient characteristics will vary between individuals; use both direct visual and audible focusing to achieve the optimum ultrasound signal.

Refer now to the appropriate Operating Handbook for positioning the Probe and operating the monitor.

The total probe usage time remaining is displayed when the Probe is connected.

**8. Removal of the Probe from the Patient:**

When temporarily withdrawing the probe from a patient (e.g., when resting or moving the patient), the probe can be disconnected from the Patient Interface Cable. If the period of removal is transitory, the probe should be cleaned using a sterile wipe before re-insertion.

**9. Probe Disposal:**

Used probes should be disposed of in accordance with the appropriate Health Authority, or similar controlling bodies, guidelines for clinical waste.

**10. Storage:**

Store in an enclosed box between -20°C and 60°C.

Do not expose to direct Ultra Violet (U.V) light, strong smells or dampness/wet.

Manufactured By Deltex Medical Ltd, Terminus Road, Chichester, PO19 8TX, UK

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**Proizvodi Deltex Medical Ltd, Terminus Road, Chichester, PO19 8TX, UK**
Consultancy Services Ltd.