(UK) - Instructions For Use (IFU): Kinder Oesophageal Doppler Probe

1. Indications – The TrueVue/ODM+ system's beat-to-beat data on cardiovascular status are used by the managing clinician to evaluate and optimise haemodynamic performance of patients who are; in surgery and are anaesthetised or sedated; in intensive care; in the emergency room; in the obstetrics department; in other wards or departments where haemodynamic measurements are required.

The paediatric nomogram is valid for patients -
Age 0 to 15 yrs.

Weight 3 kg to 60 kg (33 lb to 132 lb).
Height 50 cm to 170 cm (20" to 67").

2. Intended Use:

The purpose of the TrueVue™/ODM+™ system is to measure and calculate the haemodynamic parameters used to assess Preload, Afterload and Contractility of the heart on a beat-to-beat basis in real time. Parameters related to haemodynamics including Stroke Volume, Stroke Distance, Cardiac Output, Peak Velocity, Systemic Vascular Resistance, Mean Arterial Pressure, Pulse Pressure Variation, Stroke Volume Variation are examples of measurements and calculations to be made available to clinicians managing patient haemodynamics.

The Kinder Oesophageal Doppler Probe is for oral use in patients under general anaesthesia.

Before use read these instructions in conjunction with the Operating Handbook.

3. Clinical benefit: Proven effective for haemodynamic optimisation to reduce complications and shorten length of hospital stay

4. Warnings - Do Not Use in the Following Situations - Patient Harm May Result:

- Do not use after the expiry date printed on the shelf box and unit pack labels.
- Do not use undue force upon insertion. Remove if difficulty arises during insertion and seek advice.
- Use only by trained paediatricians or under such supervision.
- Not recommended for use in patients under going intra-aortic balloon pumping.
- The probe may be unreliable during periods of aortic cross clamping.
- Consider coagulation status of patient prior to insertion of probe.
- Down's syndrome children may have intrinsic narrowing of their hypopharyngeal structures preventing insertion of the probe.
- It is not recommended that the KDPn probes are left in place for extended periods of time.
- Any serious incident that occurs while using this product shall be reported to Deltex Medical Ltd and your member state competent authority.
- MR unsafe, do not use this equipment in the MRI scan room
- Single patient use only. Risk of cross contamination if probe is used on more than one patient

5. Cautions:

- Do not use if the pouch has been unintentionally opened, or if the probe or pouch show signs of damage.
- The probe shall be disposed of according to local procedures
- Exercise caution when handling. The probe contains an internal spring that will cause the probe to uncoil when released from the Pouch.
- Open the pouch at the chevron seal end.
- Do not use scissors or sharp tools to open the pouch.
- Clinicians are advised to take into consideration the potential of increased risks of pressure on the oropharynx when used in combination with other equipment.
- Clinicians should consider the removal of the probe between periods of measurement.
- Should the probe be removed and stored for same patient reuse, refer to local guidelines for cleaning and storage.

6. Contraindications:

- Do not use on patients weighing less than 3 kg.
- Do not insert nasally.
- Do not use with carcinoma of the pharynx; larynx; oesophagus.
- Do not use with aneurysms of the thoracic aorta.
- Do not use with proximal coarctation of the aorta.
- Do not use with tissue necrosis of the oesophagus.
- Do not use in close proximity to laser surgery
- Not for use in patients with pharyngo-oesophago -gastric pathology and/or severe bleeding diatheses

7. General Information:

The Doppler Probes are supplied STERILE.
The Probe can only be used with a Deltex Medical TrueVue/CardioQP/CardioQ ODM/CardioQP ODM Monitoring System. Probe output specifications is documented within the ODM Monitor handbook.

8. Removal of the Probe from the Unit Pack:

Open the pouch at the chevron seal end and, holding the Probe and, withdraw it from the pouch. Allow the Probe to uncoil, keeping it free from contact with any object that may compromise its state of cleanliness.

9. Insertion of the Probe into the Patient:

Insertion should be carried out by an appropriately trained and qualified paediatrician or by a person under such supervision.
The Probe must be inserted orally. Note the Warnings above.
The Probe connector should be plugged into the Patient Interface Cable connector and the correct nomogram details entered.
Apply a water-based lubricant to the tip of the Probe. This aids insertion and signal acquisition.
Insert the Probe with the bevelled edge of the tip facing the back of the throat so that when in the oesophagus, the bevelled edge faces backwards and downwards towards the descending aorta.
Excessive force should not be used when inserting the Probe. Remove if undue difficulty is experienced.
Warning - Risk of patient harm if forced.
The numbered depth markers on the Probe shaft are for use as a guide for oral placement. Signal acquisition and depth of insertion will depend upon patient height. See Table 1 as a guide to expected insertion depths; use both direct visual and audible focusing to achieve the optimum ultrasound signal. Refer to the Operating Handbook for positioning the probe and operating the Monitoring System.
The total probe usage time remaining is displayed when the probe is connected.

Table 1: Paediatric Oesophageal Doppler Probe Insertion Depth

Patient Length (cm)	Doppler Depth (cm)	Marker Position (cm)
50-55	15-17	Between 15 and 20
56-60	17-20	Between 15 and 20
61-70	16-23	Between 15 and 25
71-80	17-24	Between 15 and 25
81-100	19-28	Between 15 and 30
101-120	24-28	Between 20 and 30
121-140	27-34	Between 25 and 35
>140	27-40	Between 25 and 40

This table is created from : Use of transoesophageal Doppler ultrasonography in ventilated pediatric patients : derivation of cardiac output. Tibby SM, Hatherill M, Murdoch IA. Crit Care Med. 2000 Jun;28 (6):2045-50.

10. Removal of the Probe from the Patient:

When temporarily withdrawing the probe from a patient (e.g., when resting or moving the patient), the probe can be disconnected from the Dopplink module (TrueVue System) or Patient Interface Cable (CardioQ-ODM+). If the period of removal is transitory, the probe should be cleaned using a sterile wipe before re-insertion.

11. Probe Disposal:

Used probes should be disposed of in accordance with the appropriate Health Authority, or similar controlling bodies, guidelines for clinical waste.

12. Storage:

Store in an enclosed box between -20°C and 60°C.
Do not expose to direct Ultra Violet (UV) light, strong smells or dampness/wet.

Manufactured By Deltex Medical Ltd, Terminus Road, Chichester, PO19 8TX, UK

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Proizvod Deltex Medical Ltd, Terminus Road, Chichester, PO19 8TX, UK
Consultancy Services Ltd.
medReg GmbH